07/06/2016

07/06/2016 11:46

	Additive: DMSO	0	%		
50 uL	1X_65ul 1.5ul Mg+	1ul	х	Master Mix (no primers)	7 x
Amount	Compound		Amount	Compound	
35.5 uL	ddH_2O	0.71	249 uL	ddH_2O	248.5
0 _{uL}	DMSO 100%	0	0 _{uL}	DMSO 100%	0
5 _{uL}	10x PCR buffer	0.1	35 _{uL}	10x PCR buffer	35
4 _{uL}	dNTPs (2.5mM each)	0.08	28 _{uL}	dNTPs (2.5mM each)	28
3 _{uL}	MgCl2	0.06	21 _{uL}	MgCl2	21
0.5 _{uL}	Pf, 10μM	0.01	3.5 _{uL}	Pf, 5μM	332.5
0.5 _{uL}	Pr, 10 μM	0.01	3.5 _{uL}	Pr, 5 μM	
1 uL	DNA template	0.02	7 _{uL}	DNA template	339.5
0.5 _{uL}	Polymerase	0.01	3.5 _{uL}	Polymerase	
50 uL	Total	1	350 uL	Total mm per25 ul:	48.5

original protocol with 5x enchancer

Notes:

Genomic DNA BL2 tested in 16S PCR reaction, 2 types with and with no DMSO

94 C - 5 min

94 C, 45 s: 48 C, 1 min: 72 C, 2 min: 33x

72 C, 5 min

4 C, inf

primers done from stock: 9f 400, 1492r 200: by dilution to 50

	ul stock: uM	h2o ul resulting uM	ul resulting
9f	20	400 140 50	160
1492r	40	200 120 50	160

	Prepare Mmix Taq	0	%			
1000 uL	1X_65ul 1.5ul Mg+	1ul	Х	Master Mix (no primers)	1 x	
Amount	Compound		Amount	Compound		
690 uL	ddH_2O	0.69	690 uL	ddH_2O	690	
0 uL	DMSO 100%	0	0 uL	DMSO 100%	0	
100 uL	10x PCR buffer	0.1	100 uL	10x PCR buffer	100	
80 uL	dNTPs (2.5mM each)	0.08	80 uL	dNTPs (2.5mM each)	80	
60 uL	MgC12	0.06	60 uL	MgCl2	60	
20 uL	Pf, 5μM	0.02	20 uL	Pf, 5µM	930	
20 uL	Pr, 5 μM	0.02	20 uL	Pr, 5 μM		
20 uL	DNA template	0.02	20 uL	DNA template		
10 uL	Polymerase	0.01	10 uL	Polymerase		
1000 uL	Total	1	1000 uL	Total mm per25 ul:	930	1000