

IU 1 ug of lambda DNA in 1 hour at 37 C in 50 uL

Recommended protocol for digestion

Nuclease free water	16 uL
10X Buffer B	2 uL
DNA (0.5 -1 ug/uL)	1 uL
Apal (0.5-2 uL)	1 uL

Total: 20 uL
Mix gently and spin down for a few seconds

Incubate at 37 C for 1 - 16 hours

Recommended protocol for PCR products

PCR reaction mixture (0.5 ug of DNA)	10 uL
nuclease-free water	18 uL
10X Buffer B	2 uL
Apal (1 -2 uL)	1 uL

Total: 31 uL

Mix gently and spin down for a few seconds

Incubate at 37 C for 1 - 16 hours

Our protocol (extracted from double digest)

Nuclease free water	24.5 uL
10X Buffer B	5 uL
DNA	20 uL
Apal	0.5 uL

Total: 50 uL

Our protocol for double digest

for PCR product

Nuclease free water	23.5 uL
10X Buffer B	5 uL
DNA (0.5 -1 ug/uL)	20 uL
SacI	1 uL
Apal	0.5 uL

Total: 50 uL

for Plasmid

pCM184	43.5 uL
nuclease-free water	0 uL
10X Buffer B	5 uL
SacI	1 uL
Apal	0.5 uL

Total: 50 uL

DNA 1

conc	7.3 ug/mL
length	1085 bp
Avg Mw (ss)	329593.5 g/mol
Avg Mw (ds)	659187 g/mol
n (ss)	2.21485E-11 mol/mL

DNA2

conc	72.1 ug/mL
length	6760 bp
Avg Mw (ss)	2053091 g/mol
Avg Mw (ds)	4106182 g/mol
n (ss)	3.51178E-11 mol/mL

Ratios

length	6.2304
conc (mass)	9.8767 1.5852
conc (molar)	1.5856

50 uL						
My protocol for double digest						
for PCR product			for Plasmid			
Nuclease free water	0 uL		nuclease-free water		153 uL	
10X Buffer B	10 uL		10X Buffer B		15 uL	
PCR purified DNA	80 uL		Plasmid purified DNA		20 uL	
SacI, 10U/uL	2 uL		SacI, 10U/uL		8 uL	
Apal, 10U/uL	1 uL		Apal, 10U/uL		4 uL	
Total:	93 uL		Total:		200 uL	
Apal	10U/uL	min	5 U	for	1ug	dna
SacI	10U/uL	min	5 U	for	1ug	dna
	insert mass				Plasmid mass	
	7 ng/uL	560 ng	250 ng/uL		5000 ng	
min Apal Enzyme	0.1 uL				25 uL	
min SacI Enzyme: 2x Apal	0.2 uL				50 uL	

http://www.methods.info/Methods/RNA_DNA/restr_analysis.html

DNA up to 1ug
Enzyme 1 uL

<https://www.lifetechnologies.com/order/catalog/product/ER1411>

DNA sample up to 30% of the reaction volume

First digest (1 hour 37 C)

for PCR product		for Plasmid	
Nuclease free water	5 uL	pCM184	38 uL
10X Tango	5 uL	nuclease-free water	5 uL
DNA (0.5 -1 ug/uL)	38 uL	10X Tango	5 uL
KpnI	2 uL	KpnI	2 uL
BglII	0 uL	BglII	0 uL
Total:	50 uL	Total:	50 uL

Second digest (1 hour 37 C)

for PCR product		for Plasmid	
First digest	50 uL	First digest	50 uL
10 xTango	6.25 uL	10X Tango	6.25 uL
BglII	0.5 uL	BglII	0.5 uL
	uL		uL
	uL		uL
Total:	56.8 uL	Total:	56.8 uL

Single clone

1 - pcm184, 2 pcm184&pcr1, 3 - hypothetical pcm184&pcr1&pcr2

28/09/2015

28/09/2015 15:43

EcoRI fast digest

DD BglIII&KpnI fast digest	plasmids				
First digest: 2XSacI, 1X Tango Buffer			Master mix		13 x
nuclease-free water	6.5	uL	nuclease-free water	84.5	
Green buffer	1.5	uL	Green buffer	19.5	
Plasmid purified DNA	6	uL	Plasmid purified DNA	78	
EcoRI fast digest	1	uL			
Total:	15	uL	Per tube	14	15

Clones: 1 -12 Control: pcm184

plasmid concentration 55 ug/ml

Plasmid mass 0.33 ug

