~~	104	100	
Ub.	/01	/ ZU	116

## 06/01/2016 16:06

40

	Additive: DMSO	5	%		
40 uL	1X_65ul 1.5ul Mg+	5	Х	Master Mix (no primers)	2 x
Amount	Compound		Amount	Compound	
26.5 uL	ddH <sub>2</sub> O	0.64	53 <sub>uL</sub>	ddH <sub>2</sub> O	53
2 <sub>uL</sub>	DMSO 100%	0.05	4 <sub>uL</sub>	DMSO 100%	4
4 <sub>uL</sub>	10x PCR buffer	0.1	8 <sub>uL</sub>	10x PCR buffer	8
3.2 <sub>uL</sub>	dNTPs (2.5mM each)	0.08	6.4 <sub>uL</sub>	dNTPs (2.5mM each)	6.4
2.4 <sub>uL</sub>	MgCl2	0.06	4.8 <sub>uL</sub>	MgCl2	4.8
0.5 <sub>uL</sub>	Pf, 5µM	0.02	1 <sub>uL</sub>	Pf, 5μM	76.2
0.5 <sub>uL</sub>	Pr, 5 μM	0.02	1 <sub>uL</sub>	Pr, 5 μM	
0.5 <sub>uL</sub>	DNA template	0.02	1 <sub>uL</sub>	DNA template	
0.4 <sub>uL</sub>	Polymerase	0.01	0.8 <sub>uL</sub>	Polymerase	
40 uL	Total	1	<b>80</b> uL	Total  mm per25 ul:	38.1

## original protocol with 5x enchancer

Notes:

Genomic DNA BL2 tested in 16S PCR reaction, 2 types with and with no DMSO

94 C - 5 min

94 C, 45 s: 48 C, 1 min: 72 C, 2 min: 33x

72 C, 5 min

4 C, inf

primers done from stock: 9f 400, 1492r 200: by dilution to 50

	ul stock: uM	h2o ul resulting uM	ul resulting
9f	20	400 140 50	160
1492r	40	200 120 50	160

	Additive: DMSO	0	%				
40 uL	1X_65ul 1.5ul Mg+	5	х		Master Mix (no primers)	2 x	
Amount	Compound		Amou	nt	Compound		
28.5 uL	$ddH_2O$	0.69	57 u	ıL	ddH <sub>2</sub> O	57	
0uL	DMSO 100%	0	0 u	ıL	DMSO 100%	0	
4 uL	10x PCR buffer	0.1	. 8 u	ıL	10x PCR buffer	8	
3.2 <sub>uL</sub>	dNTPs (2.5mM each)	0.08	6.4	ıL	dNTPs (2.5mM each)	6.4	
2.4 uL	MgCl2	0.06	4.8	ıL	MgCl2	4.8	
0.5 uL	Pf, 5μM	0.02	1 u	ıL	Pf, 5μM	76.2	
0.5 uL	Pr, 5 μM	0.02	1 0	ıL	Pr, 5 μM		
0.5 uL	DNA template	0.02	1 u	ıL	DNA template		
0.4 uL	Polymerase	0.01	0.8	ıL	Polymerase		
40 uL	Total	1	80 u	иL	Total  mm per25 ul:	38.1	