#PCR1	Name	vol DNA	vol 16S_f	vol 16S_r	vol 2x SparQ	vol dH2O	Total
1	E.c Stock AB d1	5	1	1	12.5	5.5	25
2	E.c 1/10 AB d1	5	1	1	12.5	5.5	25
3	E.c 1/100 AB d1	5	1	1	12.5	5.5	25
4	E.c 1/1000 AB d1	5	1	1	12.5	5.5	25
5	E.c Stock AB d7	5	1	1	12.5	5.5	25
6	E.c 1/10 AB d7	5	1	1	12.5	5.5	25
7	E.c 1/100 AB d7	5	1	1	12.5	5.5	25
8	E.c 1/1000 AB d7	5	1	1	12.5	5.5	25
9	P&C Stock AB d1	5	1	1	12.5	5.5	25
10	P&C 1/10 AB d1	5	1	1	12.5	5.5	25
11	P&C 1/200 AB d1	5	1	1	12.5	5.5	25
12	P&C 1/4000 AB d1	5	1	1	12.5	5.5	25
13	P&C Stock AB d7	5	1	1	12.5	5.5	25
14	P&C 1/10 AB d7	5	1	1	12.5	5.5	25
15	P&C 1/200 AB d7	5	1	1	12.5	5.5	25
16	P&C 1/4000 AB d7	5	1	1	12.5	5.5	25
17	E.c Stock AB d14	5	1	1	12.5	5.5	25
18	E.c 1/10 AB d14	5	1	1	12.5	5.5	25
19	E.c 1/100 AB d14	5	1	1	12.5	5.5	25
20	E.c 1/1000 AB d14	5	1	1	12.5	5.5	25
21	P&C Stock AB d14	5	1	1	12.5	5.5	25
22	P&C 1/10 AB d14	5	1	1	12.5	5.5	25
23	P&C 1/200 AB d14	5	1	1	12.5	5.5	25
24	P&C 1/4000 AB d14	5	1	1	12.5	5.5	25
25	Contr A d1 (old 17)	5	1	1	12.5	5.5	25
26	Contr A d7 (old 18)	5	1	1	12.5	5.5	25
27	Contr A d14	5	1	1	12.5	5.5	25
28	nm1-19 native	5	1	1	12.5	5.5	25
29	nm1-19 dc0	5	1	1	12.5	5.5	25
30	nm1-19 u Soil fraction	5	1	1	12.5	5.5	25
31	ss 8u	5	1	1	12.5	5.5	25
32	ss 5u	5	1	1	12.5	5.5	25
33	ss 2u	5	1	1	12.5	5.5	25
34	ss 1u	5	1	1	12.5	5.5	25
35	ss 0.6u	5	1	1	12.5	5.5	25
36	ss 0.2u	5	1	1	12.5	5.5	25
37	ss filtrate	5	1	1	12.5	5.5	25
38	nm3-5 8u	5	1	1	12.5	5.5	25
39	nm3-5 5u	5	1	1	12.5	5.5	25
40	nm3-5 2u	5	1	1	12.5	5.5	25
41	nm3-5 1u	5	1	1	12.5	5.5	25
42	nm3-5 0.6u	5	1	1	12.5	5.5	25
43	nm3-5 0.2u	5	1	1	12.5	5.5	25
44	nm3-5 filtrate	5		1	12.5	5.5	
45	ss dc2	5	1	1	12.5	5.5	25
46	ss dc3	5	1	1	12.5	5.5	25

47	nm1-19 8u	7.5	1	1	12.5	3	25
48	nm1-19 5u	7.5	1	1	12.5	3	25
49	nm1-19 2u	7.5	1	1	12.5	3	25
50	nm1-19 1u	7.5	1	1	12.5	3	25
51	nm1-19 0.6u	7.5	1	1	12.5	3	25
52	nm1-19 0.2u	7.5	1	1	12.5	3	25
53	nm1-19 filtrate	7.5	1	1	12.5	3	25
54	1-19 dc1	7.5	1	1	12.5	3	25
55	1-19 dc2	7.5	1	1	12.5	3	25
56	ss2 8u	7.5	1	1	12.5	3	25
57	ss2 5u	7.5	1	1	12.5	3	25
58	ss2 2u	7.5	1	1	12.5	3	25
59	ss2 1u	7.5	1	1	12.5	3	25
60	ss2 0.6	7.5	1	1	12.5	3	25
61	ss2 0.2	7.5	1	1	12.5	3	25
62	ss2 filtrate	7.5	1	1	12.5	3	25
63	nm3-5 native	2.5	1	1	12.5	8	25
64	nm3-5 u Soil fraction	2.5	1	1	12.5	8	25
65	ss native	2.5	1	1	12.5	8	25
66	ss dc0	2.5	1	1	12.5	8	25
67	ss dc1	2.5	1	1	12.5	8	25
68	ss u Soil fraction	2.5	1	1	12.5	8	25
69	ss2 u Soil fraction	2.5	1	1	12.5	8	25
70	- Control	2.5	1	1	12.5	8	25

DNA5 Ma	aster Mix 50 re	eactions	per tube		
	1x	х	20 uL		
dH2O	5.5	275			
forward	1	50			
reverse	1	50			
SparQ	12.5	625			
	20	1000			

	DNA7.5 Mass	ter Mix 1 reac	tions	per tube	
	1	x x		17.5 uL	
	dH2O	3	54		
	forward	1	18		
	reverse	1	18		
	SparQ	12.5	225		
16		17.5	315		

	DNA2.5 N	Master Mix	( 10 re	actions	per tube
		1x	Х		22.5 uL
	dH2O		8	80	
8	forward		1	10	
	reverse		1	10	
	SparQ	1	2.5	125	
		2	2.5	225	