27/04/2016

27/04/2016 15:16

	Additive: DMSO	0	%			
30 uL	1X_65ul 1.5ul Mg+	1ul	Х	Master Mix (no primers)	2 x	
Amount	Compound		Amount	Compound		
20.7 uL	ddH ₂ O	0.69	41.4 _{uL}	ddH_2O	41.4	
0 uL	DMSO 100%	0	$0_{ m uL}$	DMSO 100%	0	
3 _{uL}	10x PCR buffer	0.1	6 _{uL}	10x PCR buffer	6	
2.4 uL	dNTPs (2.5mM each)	0.08	4.8 uL	dNTPs (2.5mM each)	4.8	
1.8 uL	MgC12	0.06	3.6 uL	MgCl2	3.6	
0.6 _{uL}	Pf, 5μM	0.02	1.2 uL	Pf, 5μM	55.8	
0.6 uL	Pr, 5 μM	0.02	1.2 uL	Pr, 5 μM		
0.6 uL	DNA template	0.02	1.2 uL	DNA template		
0.3 uL	Polymerase	0.01	0.6 uL	Polymerase		
30 uL	Total	1	60 uL	Total mm per25 ul:	27.9	30

original protocol with 5x enchancer

Notes: Genomic DNA BL2 tested in 16S PCR reaction, 2 types with and with no DMSO

94 C - 5 min

94 C, 45 s: 48 C, 1 min: 72 C, 2 min: 33x

72 C, 5 min 4 C, inf

primers done from stock: 9f 400, 1492r 200: by dilution to 50

	ul stock: uM	h2o ul resulting uM	ul resulting
9f	20	400 140 50	160
1492r	40	200 120 50	160

	Prepare Mmix Taq	0	%		
500 uL	1X_65ul 1.5ul Mg+	1ul	Х	Master Mix (no primers)	1 x
Amount	Compound		Amount	Compound	
345 uL	ddH ₂ O	0.69	345 uL	ddH_2O	345
0uL	DMSO 100%	0	0 uL	DMSO 100%	0
50 uL	10x PCR buffer	0.1	50 uL	10x PCR buffer	50
40 uL	dNTPs (2.5mM each)	0.08	40 uL	dNTPs (2.5mM each)	40
30 uL	MgCl2	0.06	30 uL	MgC12	30
10 uL	Pf, 5μM	0.02	10 uL	Pf, 5µM	465
10 uL	Pr, 5 μM	0.02	10 uL	Pr, 5 μM	
10 uL	DNA template	0.02	10 uL	DNA template	
5 uL	Polymerase	0.01	5 uL	Polymerase	
500 uL	Total	1	500 uL	Total mm per25 ul:	465